


# Examples of Colonic Hydrotherapy Machines

# Machines vs. Gravity

## **What is the difference between working with a CH machine and with a Gravity (tank) setup?**


- Ease of use - machines are thought to be easier – though this is debatable - possibly better in a busy clinic
- Pressure regulation – height of water vs. mechanical control
- Integrated viewing tube on some models – some prefer this
- Quicker temperature changes with a machine, but with experience a tank can be managed to achieve the therapist's/client's needs
- Price – gravity system much cheaper

# Gravity system




© ARCH Manual of Good Practice and Hygiene  
(Revised June 2022)

# Gravity system pressure



# Colonic Supplies (UK)

- Legato II
  - (Galina Imrie)
- Basic but cheap(ish)
- No waste view
- ~£3200
- [Colonic Supplies link](#)



# Dotolo (USA)

- **“Toxygen”**
  - (Ray Dotolo)
- Old-established
- Waste view tube
- No longer available new, but...



# QMS Dotolo-like machine (USA)

- “Aqua Cleanse”
- Waste view tube
- ~£9000
- CS are UK agents
- [QMS Arizona link](#)





# QMS (USA)





# Herrmann (Germany)

- **High quality**
- Waste view tube
- £8000/~£11000
- Water volume metering ('Comfort' model)
- UK support
- ["Dotolo Europe" - UK supplier](#)



# Transcom (Spain)

- **Wide range of models**
- Waste view tube
- HC-1 is manual
- HC-2000 is automatic with remote control
- €5500 - €10100
- [Transcom website](http://www.transcom.es)




# Videos

- **There are some useful videos at this site:**
  - <https://dotoloeurope.com/colon-hydrotherapy-treatment-videos/>
- Particularly:
  - Herrmann demonstration (esp. how to use the disinfectant cycle)
  - Dotolo overview [old, from their heyday, for historical interest]

# Incoming water isolation

- **Check your local water regulations**
- a “double check valve” provides isolation to EN1717 Cat 3 (~£5 x2)



# Water heating/isolation

- **Herrmann box**
- Heater
- Water isolation to EN1717 Cat 5
- £3500



# End

**Thank you**